

PHILIPS

Subject: Affiniti Analysis packages and DICOM SR (v1.5 and higher)

Affiniti uses the following designations for positioning purposes:

- “70” is the upper model that offers the most options and transducer choices
- “50” is the same physical system and software but offers fewer features and/or transducers
- “G” is optimized for General Imaging, Shared Service, and Radiology use
- “C” is optimized for Cardiology use
- “W” is optimized for Women’s Healthcare use

NOTE: The analysis library is shared between EPIQ and Affiniti platforms

Regardless of configuration, *all systems use the identical analysis packages and measurement libraries*. Because the analysis packages are the same, we use a single [DICOM Conformance Statement](#) that covers both the Affiniti 50 and Affiniti 70. Not all configurations support all DICOM SR Templates. Support is per the following table:

	OB/Gyn 5000	Vascular 5100	Adult Echo 5200	Pediatric Echo 5220	Pediatric Echo Private	Small Parts Private
Affiniti 70G	X	X	X	X	X	X
Affiniti 70C		X	X	X	X	
Affiniti 70W	X					X
Affiniti 50G	X	X	X	X	X	X
Affiniti 50C		X	X	X	X	
Affiniti 50W	X					X

